

Standardization and Metrological Traceability of a Commercial anti-MAG (Myelin Associated Protein) Antibodies ELISA

V. ECKHARDT¹, I. WINTER¹, H. ILGÜ¹, J. AFONSO¹, S. KRÄUCHI¹, R. COTTI¹, G. MILJOJKOVIC¹, B. RICKEN¹, M. SCHNEIDER¹, C.B. GERHOLD¹



¹ BÜHLMANN Laboratories, Schönenbuch, Switzerland

INTRODUCTION

The BÜHLMANN anti-MAG Antibodies ELISA (MAG-ELISA) is the acknowledged gold standard *in vitro* diagnostic ELISA, to reliably quantify anti-MAG IgM antibodies in demyelinating neuropathies, a rare autoimmune-disease. There are no recognized reference materials or reference measurement procedures for anti-MAG antibodies. We guarantee measurement consistency for the end-users over time, with a transparent traceability chain. This is achieved by using an internal reference material (IRM) to produce standardized calibrators (Fig. 3). The anti-MAG antibodies ELISA is therefore ready for IVDR, the new European regulatory basis for *in vitro* diagnostic medical devices.

METHODS

An IRM is generated from pooled individual patient sera, which are positive for anti-MAG IgM antibodies. Following the protocol by Blirup-Jensen et al., 2008 (Fig. 3), the value of the IRM is assigned to a calibrator stock, which is subsequently gravimetrically diluted into calibrators (Fig. 2). The IRM traceable MAG-ELISA was compared to the current version (Fig. 4). In addition, within-laboratory precision (total), repeatability (within run) and reproducibility (total) were assessed, measuring five samples covering the measurement range. Within-laboratory precision (total) and repeatability (within-run) were assessed with 20 days * 2 runs * 2 replicates ($n_{tot} = 80$), reproducibility (total) with 3 instruments/ lots/ operators con-founded * 5 days * 1 run * 5 replicates ($n_{tot} = 75$).


RESULTS

The combined uncertainty of the calibrator at the 95% confidence interval is 31% (Fig. 3). The IRM traceable calibrators demonstrate acceptable trueness when compared to the current calibrator material: the slope in the Passing-Bablok was 0.97 and a mean bias of 8.1% was determined for the measuring range of the MAG-ELISA (Bland-Altman analysis) (Fig 4). The within-laboratory precision (total) %CV was less than 20% with a range of 5.5 - 15.9%. The repeatability (within run) %CV was less than 15% with a range of 3.2 - 11.8% (Fig. 5). The reproducibility (total) %CV was less than 25% with a range of 10.0 - 21.6% (Fig. 6).


FIGURES

Figure 1:

a) Principle of the assay: sandwich ELISA




b) Myelin Associated Glycoprotein



Pronker et al, 2016 Nature Communications Doi: 10.1038/NCOMMS13584

Figure 2: Metrological Traceability & Standardization



CLSI Guideline EP32-R Metrological Traceability and Its Implementation; A Report


Figure 3:

a) Value Transfer Protocol


	dilutions	# runs	# replicates	# days
IRM	6	3	2	4
Calibrator Material	5 - 6	3	2	4
IRM Mastercontrol	1	3	2	4

IRM = Internal Reference Material

b) Value Assignment Calibrator Stock



c) Uncertainty Calculation



$$u_{comb} \text{ of Calibrator at } 95\% \text{ CI} = 2x \sqrt{u_{IRM}^2 + u_{Value Transfer}^2 + u_{CASET Dilution}^2} = 2x \sqrt{13^2 + 9^2 + 0.41^2} = 31\%$$

The combined uncertainty of the calibrator at the 95% CI is 31%.

Blirup-Jensen et al., 2008, DOI 10.1515/CCLM.2008.289

Figure 4: Current assay vs. IVDR assay

The IRM traceable MAG-ELISA was compared to the current version. The slope in the Passing-Bablok was 0.97 and the mean bias in the Bland-Altman was 8.1%.




Figure 5: Within-laboratory precision (total) and repeatability (within run): 20 days * 2 runs * 2 replicates (n total = 80)

Sample	Mean, BTU	N	Within-Run		Between-Run		Between-Day		Within-Laboratory	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV
2	2250.7	80	266.6	11.8	198.5	8.8	130.3	5.8	357.0	15.9
3	8849.3	80	348.7	3.9	314.3	3.6	121.9	1.4	485.0	5.5
4	19683.3	80	622.1	3.2	1492.5	7.6	908.4	4.6	1854.6	9.4
5	37185.0	80	1683.7	4.5	3082.5	8.3	1465.5	3.9	3805.9	10.2

The within-laboratory precision (total) %CV was less than 20% with a range of 5.5 - 15.9%. The repeatability (within run) %CV was less than 15% with a range of 3.2 - 11.8%

Figure 6: Reproducibility

3 instruments/ lots/ operators * 5 days * 1 run * 5 replicates (n total = 75).

Sample	Mean, BTU	N	Within-Run		Between-Day		Between-Lot-Inst.-Operator		Reproducibility	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV
2	2802.0	75	180.8	6.5	516.6	18.4	260.8	9.3	606.3	21.6
3	9051.8	75	258.2	2.9	820.7	9.1	278.7	3.1	904.4	10.0
4	18240.7	75	530.9	2.9	1146.3	6.3	1474.6	8.1	1941.8	10.6
5	34713.4	75	893.5	2.6	2740.0	7.9	2023.4	5.8	3521.4	10.1

The reproducibility (total) %CV was less than 25% with a range of 10.0 - 21.6%.

CONCLUSIONS

The transparent traceability chain of the anti-MAG Antibody ELISA not only truly displays the assay's uncertainty and leads to stable results, but also confirms its gold standard status.



Correspondence: cbg@buhlmannlabs.ch
All authors are employees of BÜHLMANN Laboratories AG.